











Suitable for people with a severe to very severe disability

6100 2500 / 6100 2510

Sirocco basic & deluxe

- · Maximum load 227 kg.
- Height adjustable, electrically operated and equipped with a removable battery, emergency descent function, emergency stop and user-friendly charger.
- · Back adjustment up to 25 degrees (Sirocco deluxe).
- · Easy and safe operation of side rails with locking.
- Corrosion resistant due to aluminum, stainless steel, plastic and galvanized parts.
- · Sleek finish for optimal cleaning and disinfection.
- · 2 position side rail (normal and wide).
- Comfortable mattress, antibacterial, flame retardant and free of harmful plasticizers.
- · Water drainage outlet hose left or right.












With built-in length adjustment (Sirocco basic & deluxe)



Sirocco deluxe & Sirocco basic









Technical Specifications

Raise/Lower Electric

Manual Operation 4-button manual operation, IPX6

• Control box 24VDC, max. 250 W IPX6

Battery 24V, 2.9 Ah IPX6

• Battery charger 100-240V AC, 50-60 Hz, max. 1A,

charging station (wall)

Art.nr.	Art. description	Outside length	Outside width	Inside length	Inside width	High/ low	Max. load	Own weight in kg
6100 2500	Sirocco shower stretcher basic	186 206	95	170 190	70-75	45-95	227	75
6100 2510	Sirocco shower stretcher deluxe	186 206	95	170 190	70-75	45-95	227	75
6720 0005	Mattress 186 basic	-	-	170	70	-	-	-
6720 0015	Mattress 206 basic	-	-	190	70	-	-	-
6720 0010	Mattress 186 deluxe	-	-	170	70	-	-	-
6720 0020	Mattress 206 deluxe	-	-	190	70	-	-	-

All measurements in this matrix are in centimeters.



