

Stylex

Effortless care, elevated comfort.

6100 3070

Stylex hydraulic shower trolley

Built for demanding care environments, Stylex combines modern functionality with exceptional durability. Designed specifically for use in humid areas such as swimming pools and hydrotherapy facilities, Stylex features a robust stainless-steel frame that is entirely rust-free. Its corrosion-resistant structure ensures lasting reliability and hygiene, even under intensive daily use in wet conditions. With smooth, easy-to-clean surfaces and an elegant, functional design, Stylex represents Lopital's commitment to creating products that perform effortlessly – wherever care is provided.

More than a shower trolley, Stylex is a professional care solution that brings together safety, comfort, and style. The height-adjustable platform allows for an ergonomic working posture, while the stable frame and premium materials guarantee confidence during every transfer. Resistant to moisture and cleaning agents, and designed for longevity, Stylex offers the perfect balance between practicality and comfort – ensuring both caregiver ease and client well-being, every day.

Ergonomic Use. Height-adjustable design and smooth positioning reduce physical strain, ensuring comfortable and efficient care at any level.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.



 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com


Technical specifications

Max. load	135 kg / 297.6 lbs
Product weight	50 kg / 110.2 lbs
Length	190 cm / 74.8 inch
Width	72 cm / 28.3 inch
Height (adjustable)	52 to 95 cm / 20.5 to 37.4 inch




Suitable for people with severely reduced to very severely reduced mobility

The water-permeable covering offers a comfortable lying position and is kind to the skin



Adjustable side rails.



Adjustable backrest.

Hydraulic height adjustment. For a healthy working posture and minimal physical strain during patient care



Stable and reliable. The sturdy frame and the smoothly rolling castors with a central brake ensure control and safety

Smooth surfaces and clever drainage ensure quick cleaning and optimal hygiene

