

6500 4025 / 6500 5030

Vienna

Suitable for people with a moderate to very severe disability




Technical specifications

Product no.	Product description	Outside length	Outside width	Inside length	Inside width	Lying surface height adjustment	Max. load kg	Max. load lb	Product weight in kg	Product weight in lbs
6500 4025	Vienna dressing table hydr 190	190	82	-	70	56 - 106	150	330	59	130
6500 5030	Vienna dressing table wall-mounted	190	91	-	70	45 - 95	227	500	70	154

All measurements in this matrix are in centimeters.

Simple operation improves ease of use.



6500 5030
Vienna dressing table
Wall-mounted



6500 4025
Vienna dressing table
hydraulic