

Tango XXL

Ease and reliability, also in demanding care situations

5100 5800

Tango XXL shower toilet chair

The Tango XXL has been specially developed to support bariatric clients with dignity and ease. With a safe working load of up to 360 kg, the chair provides comfort through its wide, ergonomically shaped seat and backrest, complemented by fold-away armrests and swing-out footrests with calf supports. These features not only make sitting and standing easier for the user but also ensure a safe and efficient workflow for the caregiver.

Versatility is at the heart of the Tango XXL. Designed for use as both a shower and toilet chair, it can be positioned directly over many toilets or used with a bedpan or collection bucket. This multifunctional approach streamlines daily care routines, while its professional-grade construction ensures the chair fits seamlessly into demanding healthcare environments.

Highly stable design. The robust frame, low centre of gravity and braked castors ensure maximum stability during transfers and care moments.

Ergonomic Use. Equipped with smart handles and four braked wheels to simplify manoeuvring, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical specifications

| | |
|----------------|-------------------|
| Max. load | 360 kg / 800 lbs |
| Product weight | 31 kg / 68.3 lbs |
| Seat width | 70 cm / 27.6 inch |
| Seat depth | 54 cm / 21.3 inch |
| Seat height | 51 cm / 20.1 inch |



Suitable for people with very slightly to severely reduced mobility

Ergonomic handle
for smooth and light
operation

Independent and
foldable armrests
for improved
accessibility

Bumper to prevent
damage to the frame



Ergonomic seat for
toilet use and bedpan.

Ergonomic seating comfort with
removable seat parts via snap-lock
system for easy cleaning

Guide rails for bedpan or toilet bucket

Directional castor for safe and efficient
mobility of the shower chair

Optional

6950 1104
Plastic bedpan and
lid.



6840 5500
Adjustable head/
neck support.



Options:

- Product no. 6950 1104 Bedpan with lid plastic
- Product no. 6975 8044 Seatbelt 2-point
- Product no. 6975 8045 Seatbelt 4-point
- Product no. 6960 8026 Comfort leg support

Technical specifications

| | |
|-------------------------|--------------------|
| Total length base frame | 97 cm / 31.1 inch |
| Total width base frame | 83 cm / 32.7 inch |
| Turning radius | 115 cm / 45.3 inch |

LFT5800-2025-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl

All of the aids produced by Lopital comply with the applicable standards for Class 1 medical products (MDR (EU) 2017/745).



ISO CERTIFIED
9001 - 13485 - 14001

MADE IN
THE NETHERLANDS



Copyright © 2025 Lopital. All rights reserved.