

Vienna

Ergonomic care, efficient space, elevated dignity

6500 5030

Vienna Deluxe wall-mounted dressing table


The Vienna range (models 5030, 5035, 4025) redefines the modern care environment with its intelligent wall-mounted design and compact footprint. By anchoring the support frame to the wall, Vienna frees up essential floor space – ideal for tight washrooms or care rooms – while delivering a robust and secure platform. Its electrically height-adjustable surface can be precisely positioned to match wheelchairs or facilitate safe transfers, improving workflow and client mobility.


This series features smooth backrest adjustment (for 5030-type models) and intuitive folding mechanisms for quick clearances when not in use. With materials selected for optimal strength and ease of maintenance, Vienna delivers a high-performance and low-profile solution tailored for demanding care settings. Its versatile integration with lifts, sliding aids, and other assistive devices ensures Vienna becomes a central, dependable element in any hygiene or dressing routine.

Ergonomic Use. Height-adjustable design and smooth positioning reduce physical strain, ensuring comfortable and efficient care at any level.

Reliable & Long-Lasting. Built to last with a 60-month warranty for long-term reliability and peace of mind.

Antibacterial Coating & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical Specifications

High / low drive	Electric (Vienna Deluxe)
Manual	Four-button hand control with hanging hook, IPX6
Control box	24VDC, max. 250 W, IPX6
Battery box	24V, 2.9 Ah, IPX6
Charger	100-240V AC, 50-60 Hz, max 1A, Wall charger



70205003
Battery box



70205004
Hand control
with battery
indication



Suitable for people with severely reduced to very severely reduced mobility

6500 5035

Vienna Basic wall-mounted fixed height dressing table



6500 4025

Vienna hydraulic dressing table



Adjustable backrest
(Vienna 4025, Vienna Deluxe 5030)



Adjustable siderails

Art.no.	Art. description	Outer length	Outside width	Inner length	Inner width	High / low	Max. load	Product weight in kg
6500 4025	Vienna hydraulic dressing table	190 cm	82 cm / 32 inch	188 cm	70 cm 27.6 inch	56-106 cm 22-41.7 inch	150 kg / 330 lbs	59 kg / 130 lbs
6500 5030	Vienna Deluxe wall-mounted dressing table	190 cm	91 cm / 35 inch	188 cm	70 cm 27.6 inch	45-95 cm 17.7-37.4 inch	227 kg / 500 lbs	75 kg / 165 lbs
6500 5035	Vienna Basic wall-mounted fixed height dressing table	190 cm	100 cm / 36.6 inch	188 cm	70 cm 27.6 inch	-	150 kg / 330 lbs	73 kg / 160.9 lbs

LFT4025-2025-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl

