

6200 4015

# Amfora

Suitable for people with a moderate to very severe disability



6200 4015  
**Amfora shower trolley**  
**hydraulic height adjustment**  
with collection tray

## Technical specifications

| Product no. | Product description                                      | Outside length | Outside width | Inside length | Inside width | Lying sur-face height adjustment | Max. load kg | Max. load lbs | Product weight in kg | Product weight in lbs |
|-------------|--|----------------|---------------|---------------|--------------|----------------------------------|--------------|---------------|----------------------|-----------------------|
| 6200 4015   | Amfora shower trolley hydraulic 190 with collection tray | 190            | 82            | 180           | 70           | 50 - 100                         | 150          | 330           | 63                   | 139                   |

All measurements in this matrix are in centimeters.