

Amfora

Simple by design, powerful in practice

6200 4015

Amfora hydraulic shower trolley

The Amfora hydraulic shower trolley by Lopital is a height-adjustable, mobile stretcher designed for safe and comfortable showering of people with reduced mobility.

Its thoughtful design ensures optimal ease of use for both the client and caregiver due to the manually adjustable backrest that enhances user comfort and allows caregivers to adapt positioning for individual needs. Foldable side rails offer added protection during transfers and washing, while maintaining easy access for caregivers. In addition, for smooth and secure mobility, it is equipped with double-lock swivel castors, providing excellent manoeuvrability and stable positioning during use.

The Amfora represents thoughtful simplicity – a solid choice that brings together comfort, ease of use and confidence.

Plug & Play Design. Ready to use out of the box – intuitive, low-maintenance, and fast to implement in any care setting.

Ergonomic Use. Equipped with smart handles and four braked castors to simplify manoeuvrability, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.




Technical specifications

Max. load	150 kg / 330 lbs
Product weight	63 kg / 138.9 lbs
Length	190 cm / 74.8 inch
Width	74 cm / 29.13 inch
Height (adjustable)	50 to 100 cm / 19.7 to 41.3 inch



Suitable for people with severely reduced to very severely reduced mobility




Adjustable backrest.

Hydraulic height adjustment.
For a healthy working posture and minimal physical strain during care taking

The water-permeable covering offers a comfortable lying position and is kind to the skin

Smooth surfaces and clever drainage ensure quick cleaning and optimal hygiene



Stable and reliable. The sturdy frame and the smoothly rolling castors with a central brake ensure control and safety



Side rail lock.



Adjustable side rails.



LFT4015-2025-01 EN



For more information, please visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl



ISO CERTIFIED
9001 - 13485 - 14001

MADE IN
THE NETHERLANDS



Scan for
information