

Reflex XL

Empowering care with strength and precision

5100 5950

Reflex XL shower toilet chair

The Reflex XL builds on Lopital's trusted Reflex concept, offering enhanced space, comfort, and adjustability for a wide range of users. Its electric height and tilt functions allow smooth, controlled movements that improve working ergonomics while maintaining a safe, comfortable position for the client. The chair's thoughtful design supports efficient hygiene care and toilet routines, blending clinical functionality with everyday ease of use.

Crafted from high-quality materials and finished with precision, the Reflex XL provides a stable, easy-to-clean platform that fits seamlessly into professional care environments. It's designed to meet the challenges of daily use with lasting comfort, intuitive operation, and Lopital's signature focus on caregiver safety and client dignity.

Ergonomic Use. Equipped with ergonomic handle and four braked castors to simplify manoeuvrability, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

Assistive Tilt Function. Thanks to the negative seat tilt, the Reflex XL offers an effective stand-up support function. This enhances user independence while reducing physical strain for the caregiver.

 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical specifications

Max. load	227 kg / 500.5 lbs
Product weight	63 kg / 138.9 lbs
Seat width	57 – 67 cm / 22.4 – 26.4 inch
Seat depth	49 cm / 19.3 inch
Seat height (adjustable)	47 – 85 cm / 20.9 – 33.5 inch



Suitable for people with moderately to severely reduced mobility

5100 5950

Reflex XL shower toilet chair



Optional

6960 8026
Comfort leg support set



Options

- 6950 1104 Bedpan with lid plastic
- 6960 8026 Comfort leg support set
- 7057 0508 Castor 125 mm Directional

Technical specifications

Total length base frame	117 cm / 46.1 inch
Total width base frame	74 cm / 29 inch
Turning radius	132 cm / 51.9 inch
Tilt setting	0 to 30 degrees backward, 10 degrees forward
Electrical power	24V
Battery capacity	7Ah
Charger (Mascot 3743)	110-240V AC, 50-60Hz, max. 0.35A, cable length 1.5 m
Head/neck support depth/height	7 / 12 cm ; 2.8 / 4.7 inch
Footrest height adjustment	8 cm / 3.2 inch

LFT5950-2026-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl



All of the aids produced by Lopital comply with the applicable standards for Class I medical products (MDR (EU) 2017/745).

Copyright ©2026 Lopital. All rights reserved.