

Timo

Ergonomic, durable and stylish in any healthcare environment

5400 6020 / 6021 / 6030 / 6035

Timo shower toilet chair

The Timo shower toilet chair from Lopital is a comfortable and reliable solution for clients with mild to moderate disabilities. Thanks to its padded seat, foldable armrests, and removable snap-lock seat, the Timo offers optimal comfort and ease of use for both clients and caregivers. The chair has an antibacterial coating and is easy to clean, contributing to a hygienic care environment. With various designs, including mobile models and versions with footrests, the Timo is flexible for use in different care situations.

Compact & Manoeuvrable Design. With its lightweight and compact frame, Timo is easy to handle in smaller spaces – ideal for home care and tighter bathrooms.

Ergonomic Use. Equipped with smart handles and four braked castors to simplify manoeuvrability, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built to last with a 60-month warranty for long-term reliability and peace of mind.

Antibacterial Coating & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

Options

- Product no. 6840 5660 Armrests with front locking set
- Product no. 6840 5661 Armrests with front dual locking set
- Product no. 6950 1104 Bedpan with lid plastic
- Product no. 6975 8044 Seatbelt 2-point
- Product no. 6960 8026 Comfort leg support set*
- Product no. 7500 5693 Raising set 10 up to 30 mm*
- Product no. 7057 0508 Castor 125 mm Directional*

*Only Timo 5400 6030 and 5400 6035



Art.no.	Art. description	Seat height	Seat depth	Seat width	Max. load in kg	Product weight in kg
5400 6020	Timo 6020 shower chair	50 / 19.7 inch	45 / 17.7 inch	46	135 kg / 297 lbs	10 kg / 22 lbs
5400 6021	Timo 6021 shower toilet chair	50 / 19.7 inch	45 / 17.7 inch	46	135 kg / 297 lbs	10 kg / 22 lbs
5400 6030	Timo 6030 shower toilet chair	54 / 21.3 inch	49 / 19.3 inch	48	135 kg / 297 lbs	15 kg / 33 lbs
5400 6035	Timo 6035 shower toilet chair	54 / 21.3 inch	49 / 19.3 inch	48	135 kg / 297 lbs	15.5 kg / 34.1 lbs

*All sizes in this table are in centimetres.



Suitable for people with very slightly to severely reduced mobility



Discover our tools in action through the Interactive Product videos on our website: www.lopital.com

5400 6035

Timo 6035 shower toilet chair

Independently foldable
armrests

Universal rail for bedpan
and / or toilet bucket



All Timo's have a snap-lock
removable seat.

Ergonomic seat with
toilet opening

Smooth-running swivel castor
for safe and efficient mobility
of the chair

5400 6020

Timo 6020 shower chair



5400 6021

Timo 6021 shower toilet chair



5400 6030

Timo 6030 shower toilet chair



LFT6020-2025-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl



ISO CERTIFIED
9001 - 13485 - 14001

MADE IN
THE NETHERLANDS



Scan for
information

All of the aids produced by Lopital comply with the applicable standards for Class 1 medical products (MDR (EU) 2017/745).

Copyright © 2025 Lopital. All rights reserved.