

Tango XXL

More space and support where it matters most

5100 5850

Tango XXL shower toilet chair

The Tango XXL has been specially developed to support bariatric clients with dignity and ease. With a safe working load of up to 360 kg, the chair provides comfort through its wide, ergonomically shaped seat and backrest, complemented by fold-away armrests and swing-out footrests with calf supports. These features not only make sitting and standing easier for the user but also ensure a safe and efficient workflow for the caregiver.

Versatility is at the heart of the Tango XXL. Designed for use as both a shower and toilet chair, it can be positioned directly over many toilets or used with a bedpan or collection bucket. This multifunctional approach streamlines daily care routines, while its professional-grade construction ensures the chair fits seamlessly into demanding healthcare environments.

Highly stable design. The robust frame, low center of gravity and braked castors ensure maximum stability during transfers and care moments.

Ergonomic Use. Equipped with smart handle and four braked wheels to simplify manoeuvring, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.

 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com



Technical specifications

Max. load	360 kg / 800 lbs
Product weight	34 kg / 75 lbs
Seat width	70 – 80 cm / 27.6 – 31.5 inch
Seat depth	53 cm / 20.9 inch
Seat height	51 cm / 20.1 inch




Suitable for people with very slightly to severely reduced mobility

Height-adjustable backrest for personalized comfort and stability

Foldable armrests with push-button operation and width adjustment

Guide rails for bedpan or toilet bucket



Ergonomical seat for toilet use and bedpan.

Ergonomic seating comfort with removable seat parts via snap-lock system for easy cleaning

Bumper to prevent damage to the frame


Directional castor for safe and efficient mobility of the shower chair

Optional

6950 1104
Plastic bedpan and lid.



6840 5500
Adjustable head/neck support.



Options:

- 6840 5500 Adjustable head/neck support
- 6950 1104 Bedpan with lid plastic
- 6975 8044 Seatbelt 2-point
- 6975 8045 Seatbelt 4-point
- 6960 8026 Comfort leg support

Technical specifications

Total length base frame	115 cm / 45.3 inch
Total width base frame	85 cm / 33.5 inch
Turning radius	127 cm / 50 inch

LFT5850-2025-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl



ISO CERTIFIED
9001 - 13485 - 14001

MADE IN
THE NETHERLANDS



All of the aids produced by Lopital comply with the applicable standards for Class I medical products (MDR (EU) 2017/745).

Copyright © 2026 Lopital. All rights reserved.