

P:\SW-Ftek\Timo en Cleo\Timo 6025\

Signed	KA	Date	02-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

Copyright reserved by law

Unit of measurement: mm


# Overview parts Timo

Lopital Nederland bv

Type nr: 54006025

## Timo 54006025

Part nr.	QTY.	Art.nr.	Part
1	1	70203501	Frame
2	1	70203509	Backrest
3	4	70519928	Castor Ø 125 mm braked
4	1	70602007	Seat
5	1	70565001	Armrest
6	4	70602006	Snaploc
7	2	75005640	Cylinder head screw M8x16 Super Duplex 4529
8	4		Cap castor
9	4		Insert Ø25
10	4		Countersunk screw M6x10
11	4		Capring M4
12	4		Chipboard screw 4,2x30
13	4		Screw cap M4



P:\SW-Tek\Timo en Cleo\Timo 6025\

Signed	KA	Date	02-03-2016
Revised	01 /KA	Date	29-10-2018
Revised	02 /KA	Date	31-03-2020

## Overview parts Timo

Copyright reserved  
by law

Unit of measure-  
ment: mm

Lopital Nederland bv

Type nr: 54006025