

Flexo

Familiar comfort designed for everyday care

5100 5199

Flexo shower toilet chair

The Lopital Flexo is a compact, versatile shower and toilet chair that excels in ergonomics, flexibility and ease of use. Thanks to its smart, open design, the Flexo offers both the client and the caregiver optimal comfort in daily care and hygiene.

The chair is height-adjustable, which ensures a healthy working posture and less physical strain for the care professional. The foldable armrests and footrests make transfers easier and safer, while the handles and manoeuvrable castors ensure smooth manoeuvrability – even in smaller spaces.

Made of durable, corrosion-resistant materials with hygienic, smooth surfaces, the Flexo is easy to clean and withstands intensive daily use. The modern, functional design fits into any professional care environment.

Reliable, stylish and well thought-out – the Flexo is designed for the practice of today and tomorrow.

Hydraulic Height & Tilt Adjustment. Smooth hydraulic height control and tilt from 0° to 20° allow ergonomic working positions and safe, comfortable user positioning.

Easy & Safe Transfers. Fold-up, lockable armrests and swing-out footrests with calf support make transfers easier while maintaining stability.

Ergonomic Use. Equipped with smart handles and four braked castors to simplify manoeuvrability, ensuring high stability on various surfaces and minimizing the strain on the caregiver's joints.

Reliable & Long-Lasting. Built for daily care and backed by a 60-month warranty, offering a dependable performance.

Antibacterial & Corrosion-Resistant. The antibacterial coating with silver-ion supports hygiene and infection control, while the corrosion-resistant frame with smooth surfaces ensures durability and easy cleaning.



 Discover our tools in action through the Interactive Product videos on our website: www.lopital.com

Technical specifications

Max. load	150 kg / 330 lbs
Product weight	38 kg / 83.8 lbs
Seat width	47 cm / 18.5 inch
Seat depth	44 cm / 17.3 inch
Seat height (adjustable)	50 to 105 cm / 19.7 to 41.3 inch



Suitable for people with very slightly to severely reduced mobility

Ergonomic handle for smooth and light operation

Armrests with front closure for safety and security

High / low control by foot pump.


Guide rails for bedpan or toilet bucket

Ergonomic seating comfort

Strong and stable frame for maximum safety and support

Footrest complete

Bumper to prevent damage to the frame



Optional

6960 8026
Comfort leg support set.



6840 5661
Dual locking armrests.



6840 5500
Adjustable head/neck support.



Options

- Product no. 6830 5650 Straight armrests
- Product no. 6840 5661 Armrests with front dual locking set
- Product no. 6950 1104 Bedpan with lid plastic
- Product no. 6975 8044 Seatbelt 2-point
- Product no. 6975 8045 Seatbelt 4-point
- Product no. 6950 5624 Platilon backrest
- Product no. 6950 5625 Platilon seat U-model
- Product no. 6960 8026 Comfort leg support set
- Product no. 7500 5693 Raising set 10 up to 30 mm
- Product no. 7057 0508 Castor 125 mm Directional

Technical specifications

Total length base frame	104 cm / 40.9 inch
Total width base frame	73 cm / 28.7 inch
Turning radius	118 cm / 46.5 inch
Tilt setting	0 to 20 degrees

LFT5199-2025-01 EN



For more information also visit our website: www.lopital.com
Lopital Nederland B.V., Laarakkerweg 9, 5061 JR Oisterwijk
Tel +31 (0)13 52.39.300, E-mail info@lopital.nl, www.lopital.nl



ISO CERTIFIED
9001 - 13485 - 14001

MADE IN
THE NETHERLANDS



All of the aids produced by Lopital comply with the applicable standards for Class 1 medical products (MDR (EU) 2017/745).

Copyright ©2025 Lopital. All rights reserved.